


# ISLA

Quick Start Guide

Sonic Face & Body Contouring Ice & Heat Roller



## A Closer Look



## Charging

### BEFORE FIRST USE

ISLA is cordless and rechargeable. Before use, charge fully for 24 hours. ISLA can run up to 40 minutes on a single charge.


### CHARGING

To recharge, connect the small end of the Charging Cord into the port on the bottom of ISLA. Then connect the USB end of the charging cord into a USB port on a computer or USB wall adapter and plug into an electrical outlet. When the ISLA's battery needs charging, the Power Indicator Light will flash red. When battery power runs out, ISLA will power off. As the battery charges, the flashing red light will turn solid red, and then green when the device is fully charged and ready for use. Recharge the ISLA only as needed for 8 hours. DO NOT overcharge which can shorten battery life. Disconnect the USB Cord from the power source when not charging. DO NOT charge ISLA for more than 30 continuous hours.

### ENHANCE YOUR TREATMENT

Remove treatment globes by gently pulling without twisting. You may place these in a freezer for 60 minutes to assist in a cooling treatment, let rest for 5 minutes before use. If a warm treatment is preferred, please place in a cup of warm/hot water for 20 minutes and then place on the ISLA.


**CAUTION: DO NOT OVERHEAT THE GLOBES. NEVER PLACE ISLA'S STAINLESS STEEL GLOBES IN THE MICROWAVE OVEN!**



See Expert Tips section for more detailed info on making the most of the Detachable Globes of ISLA.

## Usage Guide

### READY TO USE



Once you have assured of the full charge of your ISLA device you can proceed with your treatment. See the Expert Advice section for more details.

1. Cleanse skin on the face and neck, rinse clean and pat dry. You are ready to begin your treatment.

### STARTING THE TREATMENT

2. Turn on your ISLA by pressing power button on the control handle. You may select low or high speed depending on personal preference and area treated. We recommend the lower speed around the eyes and where thinner skin is located.



3. Position the ISLA onto the area of clean, dry skin. Starting from chin on your jawline, roll upwards and outwards toward your ear. Repeat up to 5 rolls in each area. Then roll up and out under your cheek bone then along it toward your ear. Roll from your nose toward your ear for the lip area. Finally, roll flat along the brows, then upwards on the forehead. Neck area can be massaged in an upward movement if desired as well as other areas of the body. **For best results, DO NOT** roll the device up then down or back and forth. Always use an upward and outward movement for best contour lifting results and to help promote lymphatic drainage of excess toxins.



## Features and Benefits

### SONIC TECHNOLOGY


ISLA is equipped with Sonic Technology and two Medical-Grade Stainless-Steel Treatment Globes for Face and Body Treatments and designed to replicate the spiral kneading action of an Esthetician's fingers.

The sonic technology has more than 30,000 movements per minute. It gently massages, tones your skin and increases circulation.

ISLA works through a unique combination of sonic vibration and manual acupressure by rolling across the face. Note, the sonic vibration may cause the rollers to rotate slightly on their own. This is normal and not a feature of the treatment.

### DETACHABLE STAINLESS-STEEL TREATMENT GLOBES

ISLA's Detachable Stainless Steel Treatment Globes accommodate different contours of the face and body providing comfortable and effortless treatment to all areas. This allows for stimulating massage at different angles to further increase the lifting and firming effects.



### EASY CLEANING OF THE TREATMENT GLOBES

Detachable Stainless-Steel Treatment Globes make cleaning a breeze! When cleaning the Globes, simply spray with 70% rubbing alcohol and wipe clean. DO NOT submerge the device in water.

## Expert Advice

ISLA's Treatment Globes can provide heat therapy and cold therapy to enhance any facial massage. The cold therapy helps to firm, calm and cool your skin, while also de-puffing your under eyes. The heat therapy can relieve aches and facial muscle tension and relax. Which help to reduce the appearance of fine lines and wrinkles.

ISLA is an ideal anti-aging device that pairs perfectly with your favorite serum, creams, moisturizers, oils & sheet masks. Using ISLA to apply your skin care products helps to ensure maximum absorption. Try using heat therapy while applying. You will get the most from your products while using the heat therapy combined with the ISLA. Apply a few drops of serum, creams, moisturizers and oils or roll over sheet mask to your clean, dry face, then gently glide Treatment Globes across your skin until the product is absorbed.

For best results, DO NOT roll the device up then down or back and forth. Always use an upward and outward movement for best contour lifting results and to help promote lymphatic drainage of excess toxins.

You may safely use the device any time of day as often as you would like. Use with your serums and/or creams in the a.m. and p.m., with heat therapy or cold therapy. The possibilities are endless in achieving firm, smooth and rejuvenated skin with ISLA!

## Additional Info

The Operating Manual in its entirety including safety, care and warranty information is available at [www.islamanual.com](http://www.islamanual.com)

Email customer service at [customerservice@spasciences.com](mailto:customerservice@spasciences.com)

© 2020 Spa Sciences. All rights reserved.

## Warnings & Warranty

### Warnings:

1. ISLA is to be used on uncompromised treatment areas. Avoid use in treatment areas where any of the following is present; any area that was recently treated with chemical peels or lasers. Avoid treatment areas where cancer is present or suspected to be present. **DO NOT** use over compromised skin. Compromised skin may include but is not limited to sunburns, areas of healing, blisters, Rosacea, Psoriasis, abrasions, cuts and acne.
2. If you have any medical concerns about using ISLA, please consult your dermatologist before use.
3. ISLA is not to be used on children. Children should not play with the device.
4. If any of the components are damaged (ISLA Device Handle, or USB Charging Cord) **DO NOT** use. **DO NOT** attempt to open or services components. Email Spa Sciences customer services at [customerservice@spasciences.com](mailto:customerservice@spasciences.com).
5. ISLA is intended only for use described in the Operating Manual. Email Spa Sciences customer services at [customerservice@spasciences.com](mailto:customerservice@spasciences.com).
6. **DO NOT** use or store ISLA anywhere that it can fall or be dropped into water. If the ISLA Device Handle, USB Charging Cord falls into water, **DO NOT** reach for it as electrocution may occur. **IMMEDIATELY**, unplug from the wall. Before charging always check to make sure all components are free from water.
7. Only use the supplied USB Charging Cord to charge your ISLA Device Handle.
8. **DO NOT** overcharge which can shorten battery life. Disconnect the USB Cord from the power source when not charging the device.

The complete Operating Manual is available online at [www.islamanual.com](http://www.islamanual.com).

For home use. Not intended for use in medical or industrial applications.

This product has been tested and found to comply with the limits for a Class B digital device pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This product generates, uses and can radiate radio frequency (RF) energy and, if not installed and used in accordance with the instructions may cause harmful interference to radio communications devices (such as televisions, baby monitors, cordless phones, etc). However, there is no guarantee that RF interference will not occur in a particular installation. If this product does cause harmful interference to radio or television reception, which can be determined by turning the product off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antennae of the other radio communications device
- Move the equipment farther away from the radio or tv
- Connect the product into an outlet on a circuit different from that to which the other device is connected
- Consult the manufacturer's customer service or an experienced radio/tv technician for more suggestions.

### Warranty:

**Your warranty is not effective unless you register the Product within 90 days of purchase. To register go online to [www.warrantyspasciences.com](http://www.warrantyspasciences.com) and complete the information required for registration.**

For a period of two (2) years from the date of purchase, Spa Sciences (the "Company") warrants the ISLA Advanced Skin care Delivery System (the "Product") against damage, defects or failure resulting from ordinary consumer use under normal home conditions according to the instructions included with the Product (a "Covered Defect"). This warranty DOES NOT include damage, defects or failure caused by misuse or improper care of the Product; impact or insertion of foreign objects into or on the Product; application of harmful chemicals or abrasives to the Product; failure to follow manufacturer's instructions; normal wear and tear; alteration or overuse of the Product; or any professional or commercial use of the Product. If there is a Covered Defect during the two (2) years warranty period the Company agrees to exchange the Product with a new equivalent Product. This is the sole remedy available for breach of this limited warranty. Any claim for breach of this limited warranty must be commenced within two (2) years following the original purchase. This limited warranty applies only to a purchase of the Product from or an authorized reseller, is valid for the original end user of the Product only and cannot be transferred under any circumstances.

THE PROVISIONS OF THE FOREGOING WARRANTIES ARE IN LIEU OF ANY OTHER WARRANTY, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL (INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE). IN NO EVENT SHALL THE COMPANY BE LIABLE FOR ANY SPECIAL, INCIDENTAL, CONSEQUENTIAL, PUNITIVE, OR INDIRECT DAMAGES BASED UPON ANY LEGAL THEORY. SOME STATES DO NOT ALLOW LIMITS ON WARRANTIES, OR ON REMEDIES FOR BREACH IN CERTAIN TRANSACTIONS. IF APPLICABLE LAW DOES NOT PERMIT THE DISCLAIMER OF IMPLIED WARRANTIES, THEN THE DURATION OF SUCH IMPLIED WARRANTIES, ARE EXPRESSLY LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY SET FORTH HEREIN.